

Your practice is *our* inspiration.®

AlgiNot | The alginate alternative



AlgiNot

A time-saving, cost-effective alternative to traditional alginate



Kerr continues to deliver an alginate alternative with outstanding value and clinical performance—now in a faster set time to increase patient comfort without sacrificing working time, plus an improved thixotropic consistency. AlgiNot™ is ideal for preliminary impressions, provisional impressions, orthodontic models, bite registration, sports guards, bleaching trays and more.


*Inspired by:
Faster, easier preliminary impressions for the dental assistant was the impulse. AlgiNot FS reduces time constraints associated with regular alginate and eliminates the mess.*

Research

As traditional alginates benefit from inherent hydrophilicity, alginate alternatives must behave the same way. AlgiNot breaks the surface tension of intraoral fluids and displaces them to register outstanding detail. And with superior elasticity, AlgiNot rebounds from deformation far better than traditional alginates—or even other alternatives—to produce an impression that is a true negative replication of the dentition.



Internal data. Available upon request.




Features

Improved thixotropic consistency. Non-slumping. Does not run. Increases patient comfort.

High deformation recovery. Does not distort upon removal or after pouring multiple casts.

Low dimensional change. Retains its shape upon set. No special storage required.

Can be disinfected. Eliminates cross-contamination concerns without affecting dimensional stability.



Internal data. Available upon request.

Ordering

AlgiNot

33034 Cartridge Intro Kit

Contains: 6 Cartridges (50 ml each), 12 mixing tips, 1 Technique Card.

33035 Cartridge Refill

Contains: 24 Cartridges (50 ml each), 1 Technique Card.

34681 Volume Intro Kit

Contains: 1 Cartridge (380 ml), 10 Mixing Tips, 1 Bayonet Ring, 1 Technique Card.

34680 Volume Refill

Contains: 6 Cartridges (380 ml each), 3 Bayonet Rings, 1 Technique Card.

AlgiNot FS

33816 Cartridge Intro Kit

Contains: 6 Cartridges (50 ml each), 6 mixing tips, 1 Technique Card.

33817 Cartridge Refill

Contains: 24 Cartridges (50 ml each), 1 Technique Card.

34683 Volume Intro Kit

Contains: 1 Cartridge (380 ml), 10 Mixing Tips, 1 Bayonet Ring, 1 Technique Card.

34682 Volume Refill

Contains: 6 Cartridges (380 ml each), 3 Bayonet Rings, 1 Technique Card.

Accessories

34746 Dynamic Mixing Tips for Volume (50 tips, 2 bayonet rings)

29333 Mix Tips, Large (pkg of 48)

28747 Extruder Gun

33171 Volume Mixer



©2009 Kerr Corporation. Lit. No. 34087A
All trademarks are property of their respective owners.